

Humaneway autoclavable rodent restrainers for mice and rats, six sizes covering 15 g to 520 g animals, with a bell-shaped entry, detachable tail slot, adjustable head stopper, and a heavy stainless-steel base for stable tail vein procedures.

SPECIFICATIONS

Mouse sizes (ID × length)	28 mm × 110 mm 32 mm × 110 mm 35 mm × 126 mm
Mouse weight range	15–30 g 20–40 g 30–50 g
Rat sizes (ID × length)	47 mm × 176 mm 64 mm × 215 mm 70 mm × 215 mm
Rat weight range	100–160 g 160–450 g 400–520 g
Base material	Stainless steel (bottom-heavy)
Sterilisation	Fully autoclavable (standard steam)

FEATURES

- ▶ Six sizes: 3 for mice (28, 32, 35 mm ID) and 3 for rats (47, 64, 70 mm ID)
- ▶ Bell-shaped tube entry for smooth, injury-free animal loading
- ▶ Detachable tail slot for unobstructed tail access
- ▶ Adjustable head stopper accommodates varying body lengths
- ▶ One-way entry/exit design minimises escape attempts
- ▶ Bottom-heavy stainless-steel base for maximum bench stability
- ▶ Silicone grip for secure handling during procedures
- ▶ Fully autoclavable, suitable for GLP environments
- ▶ Designed for use with Humaneway Tail Vein Simulators and Laboratory Needles

APPLICATIONS

- ▶ Tail vein IV injection in mice and rats
- ▶ Tail vein blood collection
- ▶ Intravenous drug dosing in preclinical studies
- ▶ Any procedure requiring tail access or temporary rodent immobilisation
- ▶ Training use with Humaneway Tail Vein Simulators